


# High sensitive robust and easy-to-use transporter assays



## Influx and efflux transporter assays




**Products** Vesicular transporter assays with efflux transporters of interest.  
 Transient over-expression in HEK293 cells

Vesicles: shape and size

Vesicles diameter  
 ~ 200 nm (DLS)  
 (PDI ~ 0.27)



Vesicular uptake assay



### Ready-to-use vesicular transporter assays

MDR1/P-gp	BCRP	BSEP
MRP1	MRP2	MRP3
MRP4	MRP5	MRP8

### Contact

Errol Byron  
 Kloosterstraat 9, 5349 AB,  
 Oss Pivot Park, the Netherlands  
 E: errol.byron@cell4pharma.com  
 www.cell4pharma.com

## Key publications

Lempers VJ, *et al.* An`microb Agents Chemother. 2016.

Gozalpour E, *et al.* Toxicol In Vitro. 32:138-45, 2016.

Te Brake LH, *et al.* Tuberculosis. 96:150-7, 2016.

Gozalpour E, *et al.* Eur J Pharmacol. 5;744:18-27, 2014.

Gozalpour E, *et al.* Mol Pharm. 2;11(6):1844-55, 2014.

El-Sheikh AA, *et al.* Transl Res. 164(1):46-56, 2014.

van Beusekom CD, *et al.* BMC Vet Res. 20;9:259, 2013.

Gozalpour E, *et al.* Toxicol Sci. 131(2):502-11, 2013.

Wiegen HG, *et al.* Mol Pharm. 7;9(5):1351-60, 2012.

Wiegen HG, *et al.* Drug Metab Dispos. 39(7):1294-302, 2011.